

Tubon

Tubon Tubular Gauze Bandage



Tubon is used in the following applications:

- Casting - Applied as the primary layer on limbs, prior to the application of cast padding and casting materials.
- Oedema Control - Applied as the primary layer on limbs prior to the application of pressure bandaging.
- Physiotherapy & Rehabilitation - The lining of splints, braces and orthoses.
- Woundcare & Dermatology - Applied over dressings and wet wraps, to prevent clothing and bedding from being soiled by ointments and creams.

Tubon offers the following benefits:

TUBON is made from Pure Cotton. Cotton is the most comfortable fibre, as it absorbs perspiration away from the skin.

TUBON allows good air circulation to the skin, due to the open knit structure.

TUBON is available in 8 sizes, which will cover fingers, limbs and torsos.

TUBON is Australian Made, by an Australian owned Company

The smaller sizes are used for dressing retention on fingers and toes

The larger sizes are used on limbs and torso



Finger Applicator

Code	Size	Width mm	Application
6000	0	14	Small Toes & Fingers
6001	1	18	Toes & Fingers
6002	2	25	Finger & Toes Over Dressings
6003	3	40	Child Arms, Hands & Feet
6004	4	65	Adult Arms, Small Legs & Hands
6005	5	85	Adult Legs, Small Thighs & Head
6006	6	110	Large Thighs, Small Trunks
6007	7	155	Medium Trunks
6008	8	240	Adult Trunks
APPL-1	Pkt 10	Plastic Applicator for Finger Size	