



Ortho Thumb CMC Brace

CMC 1 Joint Saver



The Push Ortho CMC s brace stabilises the basal joint of the thumb (CMC-1) and places the thumb in a functional position.

The functional position helps guarantee optimal hand function and reduces such symptoms as pain and loss of strength.

Features

- Seamless Brace with rounded edges.
- Simple fastening system for easy application and removal
- Lightweight but durable construction.
- Antimicrobial, can be used in water, and dried easily.

Indications:

Osteoarthritis of the CMC 1 joint
Post Operative after Care
Instability of the CMC.
Gamekeepers /Skiers Thumb

Push Ortho CMC Joint Saver Brace			
Size	Left	Right	Wrist Circumference
1	PCMC-LS	PCMC-RS	16 -19.5cm
2	PCMC-LM	PCMC-RM	19.5 - 22.5cm
3	PCMC-LL	PCMC-RL	22.5 - 26cm

